

## PYCNOMETER with THERMOMETER

Code : W013



- A pycnometer with a thermometer in engineering testing laboratories is used to measure the density or specific gravity of liquids, powders, or solids. The built-in thermometer allows for precise temperature control, ensuring accurate density measurements, as density can vary with temperature. This is essential in material characterization and quality control.

### ORDERING INFORMATION

Item	Code
PYCNOMETER with THERMOMETER - 25 ml	W013G025H
PYCNOMETER with THERMOMETER - 50 ml	W013G050H
PYCNOMETER with THERMOMETER - 100 ml	W013G100H