

LAB TOOLS > GLASSWARE

## **PORCELAIN MORTAR and PESTLE**

Code: W026



 A porcelain mortar and pestle in engineering testing laboratories are used for grinding, crushing, and mixing solids into fine powders. The non-porous nature of porcelain ensures that no material is absorbed, making it ideal for preparing samples for chemical or material analysis without contamination.

## **ORDERING INFORMATION**

Item	Code
PORCELAIN MORTAR and PESTLE - Ø 80 mm	W026C080H
PORCELAIN MORTAR and PESTLE - Ø 100 mm	W026C100H
PORCELAIN MORTAR and PESTLE - Ø 130 mm	W026C130H
PORCELAIN MORTAR and PESTLE - Ø 160 mm	W026C160H



## **W026 DATASHEET**

PORCELAIN MORTAR and PESTLE - Ø 190 mm	W026C190H
PORCELAIN MORTAR and RUBBER-HEADED PESTLE	W027P001H
GLASS MORTAR and PESTLE - Ø 100 mm	W028G100H
GLASS MORTAR and PESTLE - Ø 120 mm	W028G120H
GLASS MORTAR and PESTLE - Ø 150 mm	W028G150H
GLASS MORTAR and PESTLE - Ø 180 mm	W028G180H

## **OTHER PHOTOS**

