


JAW CRUSHER

Code : M412



- Used to crush samples with medium-hard, hard, brittle and tough materials. Efficiency and safety makes it ideal for sample preparation in laboratories and industrial plants.
- Application examples are so wide such as construction materials, glass, granite, minerals, quartz, rocks, silicon, slag, coal, oxide ceramics ... etc.
- The materials after crushing have high degree of size reduction and high fineness.
- Double action of jaws that opens/closes from the bottom and the crushing mechanism is from top to bottom.
- Supplied complete with electrical motor.

TECHNICAL SPECIFICATIONS

- Double Action of jaws:
 - Opening /closing at bottom

- Crushing from up to bottom
- Jaw opening at the bottom is adjustable.
- 380 V / 50 Hz

ORDERING INFORMATION

Item	Code
JAW CRUSHER (Small)	M412X001U
JAW CRUSHER (Large)	M412X002V
JAW CRUSHER (Extra Large)	M412X003V
SPARE JAWS for M412	M412P001H
SPARE SIDE PLATES for M412	M412P002H

OTHER PHOTOS

