

## VIBRATORY DISC MILL

Code : M430



- Used for quick grinding, assuring that no sample loss can be experienced.
- It has a very wide range grinding and finishes the job in an extremely short time.
- The maximum materials feed size is 15 mm and grinds to a final fineness of 50 micron.
- The jar size that is supplied with the equipment is 100 ml.

### **TECHNICAL SPECIFICATIONS**

- Suitable for quick grinding
- No sample loss
- Wide range grinding
- Extremely short grinding time
- Material feed size < 15 mm
- Final Fineness < 50 mic

- Jar Size 100 ml
- 380 V / 50 Hz

**ORDERING INFORMATION**

| Item Name                                      | Item Code |
|------------------------------------------------|-----------|
| JAR and SPARE DISCS for M-430                  | M430P001H |
| BODY and LOCKING MECHANISM for M-430           | M430P002H |
| SPARE MOTOR for M430                           | M430P003H |
| TUNGSTEN CARBIDE JAR and SPARE DISCS for M-430 | M430P004H |
| VIBRATORY DISC MILL                            | M430X00XU |