


JAW CRUSHER

Code : M412



- Used to crush samples with medium-hard, hard, brittle and tough materials. Efficiency and safety makes it ideal for sample preparation in laboratories and industrial plants.
- Application examples are so wide such as construction materials, glass, granite, minerals, quartz, rocks, silicon, slag, coal, oxide ceramics ... etc.
- The materials after crushing have high degree of size reduction and high fineness.
- Double action of jaws that opens/closes from the bottom and the crushing mechanism is from top to bottom.
- Supplied complete with electrical motor.

TECHNICAL SPECIFICATIONS

- Double Action of jaws:
 - Opening /closing at bottom

- Crushing from up to bottom
- Jaw opening at the bottom is adjustable.
- 380 V / 50 Hz

ORDERING INFORMATION

Item Name	Item Code
SPARE JAWS for M412X001U	M412P001H
SPARE SIDE PLATES for M412X001U	M412P002H
SPARE JAWS for M412X001U	M412P003H
JAW CRUSHER (100 Series)	M412X001U
JAW CRUSHER (170 Series)	M412X002V
JAW CRUSHER (220 Series)	M412X003V
JAW CRUSHER (500 Series)	M412X004V
JAW CRUSHER (300 Series)	M412X005V

OTHER PHOTOS

