

CONSTANT HEAD PERMEABILITY SET

Code : T252



- Used to determine the permeability of granular-grained soils

STANDARDS

ASTM D2434 • BS 1377-5 • AASHTO T215

TECHNICAL SPECIFICATIONS

Model	T252X075H	T252X114H
Permeameter stand with 4 tubes		
Transparent Plastic Permeability Cell	Ø 75 mm	Ø 114 mm
Water Reservoir Tank		
Take-Off Points	3	12 (6 blacked off)

ORDERING INFORMATION

Item	Code
CONSTANT HEAD PERMEABILITY SET - Ø 75 mm	T252X075H
CONSTANT HEAD PERMEABILITY SET - Ø 114 mm	T252X114H
MANOMETER STAND with 4 TUBES (Constant Head)	T252P001H
RESERVOIR TANK (Constant Head)	T252P002H
MANOMETER TUBE - Ø 6 mm	T252P003H
MANOMETER TUBE - Ø 8 mm	T252P004H
MANOMETER TUBE - Ø 14 mm	T252P005H
MANOMETER TUBE - Ø 22 mm	T252P006H
PERMEABILITY CELL (Constant Head) - Ø 75 mm	T252P007H
PERMEABILITY CELL (Constant Head) - Ø 114 mm	T252P008H