

ASPHALT > SURFACE PROPERTIES

ROLLING STRAIGHTEDGE, TRAVELLING BEAM (High-Low Detector)

Code : A069



- Used for checking planeness/irregularities in concrete and bituminous road surfaces.
- The travelling beam is 3 meters long.
- Deviation of the surface is shown on a scale;
 - in increments of 2 mm from 0 - 10 mm
 - in increments of 5 mm from 10 - 25 mm

STANDARDS

British Department of Transport

TECHNICAL SPECIFICATIONS

- Length : 3 m
- Equipped with scale (up to 25 mm)

ORDERING INFORMATION

Item	Code
ROLLING STRAIGHTEDGE, TRAVELLING BEAM (High-Low Detector) - 3 m	A069X03XH