

ASPHALT > BITUMEN PROPERTIES

## FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)

Code : A470



- Used to determine the maximum temperature bitumen can be safely heated without the danger of an instantaneous flash in the presence of an open flame.
- Supplied with an electronic heater, a thermoregulator, two thermometers, and a sample cup.

### **STANDARDS**

ASTM D117 • ASTM D92 • AASHTO T-48

### **SUPPLIED WITH**

- Electric heater

- Thermoregulator
- Two Thermometers
- Sample Cup

## **ORDERING INFORMATION**

| Item                                                          | Code      |
|---------------------------------------------------------------|-----------|
| FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)           | A470H00XU |
| AUTOMATIC FLASH and FIRE POINT APPARATUS (Cleveland Open Cup) | A470A00XU |