

ASPHALT > BITUMEN PROPERTIES

## FLASH and FIRE POINT APPARATUS (Tag Open Cup)

Code : A471



Used to determine open-cup flash points of volatile flammable materials having flash points between 0 and 175 °F.

### **STANDARDS**

ASTM D1310 • ASTM D3143

### **ORDERING INFORMATION**

| Item                                          | Code      |
|-----------------------------------------------|-----------|
| FLASH and FIRE POINT APPARATUS (Tag Open Cup) | A471H00XU |