

CORE CUTTER

Code : K063



- Used to cut core samples.
- Equipped with a disc holder for a maximum diameter of 230 mm.
- The distance of the disc is adjustable according to the user's needs.
- Supplied with a water pump for cooling the materials while cutting.

TECHNICAL SPECIFICATIONS

- Used for cutting core samples.
- Disc diameter: Maximum 230 mm.
- Supplied with water pump for cooling.
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

Item	Code
CORE CUTTER	K063X00XH
CORE SAMPLE MEASUREMENT HOLDER	K063P00XH