

MULTI-DECK GRADING MACHINE

Code : M042



- Used specially in mining labs where large quantities of sieving is required.
- Two sieves with decks, each of 300 x 500 mm size, and the openings are 10 mm and 5 mm unless otherwise specified during the order.
- Plastic boxes, as collectors, are supplied for each deck.

TECHNICAL SPECIFICATIONS

- Equipped with pan at the bottom.
- Adjustable sieving angle.
- 380 V / 50 Hz

ORDERING INFORMATION

Item	Code
MULTI-DECK GRADING MACHINE	M042X00XU