


SOIL > CALIFORNIA BEARING RATIO

FIELD CBR SET

Code : T002



- The California Bearing Ratio (CBR) is a penetration test carried out to evaluate the mechanical strength of road subgrades and base-coarse layers.
- Equipped with a mechanical jack (2-speed), and a load ring.

STANDARDS

ASTM 1883 • BS 1377-4 • EN 13286-47

SET

- Load Ring Capacity : 50 kN
- Mechanical Jack - 2 Speed
- CBR Penetration Piston

- Penetration Dial (2 ea)
- Dial Holder
- Extension Rods Surcharge Weight

ORDERING INFORMATION

Item	Code
LAB. CBR APPARATUS - TFT	T001A05XU
LAB. CBR APPARATUS - Load Ring	T001L05XU
LAB. CBR APPARATUS - Digital	T001M05XU
SPARE LOAD/DEFORMATION INDICATOR for T001M	T001P001X
CBR PENETRATION PISTON	T001P002X
FIELD CBR SET	T002L05XH
CONVERSION FRAME for T-001/F	T002P001H
LAB. CBR APPARATUS - Digital - Multi-Speed	T003M05XU