

SOIL > ATTERBERG LIMITS

PLASTIC LIMIT TEST SET

Code : T035



- Used to determine the plastic limit, which is the minimum moisture content at which a sample can be rolled into threads with a 3 mm diameter without experiencing breakage or cracking neither longitudinally or transversely

STANDARDS

ASTM D4318 • AASHTO T90 • BS 1377:2 • NF P94-051

SET

Item

STAINLESS STEEL ROD (Ø 3 mm)

GLASS PLATE - 300 x 300 mm

MOISTURE TIN, with LID - Small - Ø 55 x 35 mm

EVAPORATING DISH [300 ml - Ø 125 mm]

SPATULA - 100 mm

BAG for T035

ORDERING INFORMATION

Item	Code
PLASTIC LIMIT SET	T 035 X00X H
BAG for T-035	T 035 P001 H
STAINLESS STEEL ROD (Ø 3 mm)	T 035 P002 H
GLASS PLATE - 300 x 300 mm	W 153 G300 H
MOISTURE TIN, with LID - Small - Ø 55 x 35 mm	G 720 X001 H
EVAPORATING DISH [300 ml - Ø 125 mm]	W 137 C300 H
SPATULA - 100 mm	G 717 X100 H