

# DISTILLING FLASK

Code : W004



- A distilling flask in engineering testing laboratories is used for heating liquids to separate components based on their boiling points through distillation. It is commonly used in chemical analysis, purification processes, and the production of solvents or distilled water in the lab.

## ORDERING INFORMATION

| Item                       | Code      |
|----------------------------|-----------|
| DISTILLING FLASK - 100 ml  | W004G100H |
| DISTILLING FLASK - 250 ml  | W004G250H |
| DISTILLING FLASK - 500 ml  | W004G500H |
| DISTILLING FLASK - 1000 ml | W004G01KH |

---

**OTHER PHOTOS**

