

PYCNOMETER (Gay Lussac)

Code : W014



- A pycnometer (Gay Lussac type) in engineering testing laboratories is used to determine the specific gravity or density of liquids or powders. It consists of a precisely calibrated glass container, often with a stopper, ensuring accurate measurements by accounting for temperature variations during testing. This is essential for density and quality control in chemical, material, and environmental testing.

ORDERING INFORMATION

Item	Code
PYCNOMETER (Gay Lussac) - 5 ml	W014G005H
PYCNOMETER (Gay Lussac) - 10 ml	W014G010H
PYCNOMETER (Gay Lussac) - 25 ml	W014G025H

PYCNOMETER (Gay Lussac) - 50 ml W014G050H

PYCNOMETER (Gay Lussac) - 100 ml W014G100H

PYCNOMETER (Gay Lussac) - 250 ml W014G250H

PYCNOMETER (Gay Lussac) - 500 ml W014G500H

PYCNOMETER (Gay Lussac) - 1000 ml W014G01KH

PYCNOMETER (Gay Lussac) - 2000 ml W014G02KH

OTHER PHOTOS

