

PLASTIC MEASURING CYLINDER

Code : W032



- A plastic measuring cylinder in engineering testing laboratories is used for accurately measuring the volume of liquids. It is lightweight and durable, making it suitable for handling a variety of chemicals or materials, especially when glass containers might be too fragile or heavy. It provides precise measurements for experiments and sample preparation.

ORDERING INFORMATION

Item	Code
PLASTIC MEASURING CYLINDER - 10 ml	W032P010H
PLASTIC MEASURING CYLINDER - 25 ml	W032P025H
PLASTIC MEASURING CYLINDER - 50 ml	W032P050H
PLASTIC MEASURING CYLINDER - 100 ml	W032P100H

PLASTIC MEASURING CYLINDER - 250 ml

W032P250H

PLASTIC MEASURING CYLINDER - 500 ml

W032P500H

PLASTIC MEASURING CYLINDER - 1000 ml

W032P01KH

PLASTIC MEASURING CYLINDER - 2000 ml

W032P02KH