

CRADLE

Code : G100P001H



- The Cradle is designed to suspend specimens beneath a weighing balance within a Specific Gravity Frame, facilitating accurate mass measurements in both air and water for specific gravity calculations.
- Constructed from high-quality stainless steel rod, ensuring durability and resistance to corrosion during frequent immersion processes.
- Compatible with ALFA's Specific Gravity Frame, seamlessly integrating into existing laboratory setups for aggregate testing.

STANDARDS

ASTM C127 • AASHTO T85 • BS 812:2 • EN 12390-7 • EN 1097-6

ORDERING INFORMATION

Item	Code
CRADLE	G100P001H