

## SPECIFIC GRAVITY PYCNOMETER

Code : G103



- The Pycnometer is designed for accurate determination of the specific gravity and water absorption of fine aggregates, essential for quality control in construction materials.
- Consists of a durable glass jar with a non-corrodible conical brass top featuring a 10 mm hole to facilitate proper filling and venting during testing.
- Includes a rubber gasket to ensure an airtight seal, preventing water or fine particles from leaking during the testing process.
- Complies with ASTM C128 and AASHTO T84 standards, ensuring reliable and repeatable results in line with internationally recognized testing methods.

### **STANDARDS**

ASTM C128 • AASHTO T84

**TECHNICAL SPECIFICATIONS**

- Material: Glass jar with brass top
- Top Opening Diameter: 10 mm
- Capacity: 1 lt

**EQUIPPED WITH**

- Rubber gasket for airtight seal

**ORDERING INFORMATION**

Item	Code
SPECIFIC GRAVITY PYCNOMETER	G103X00XH