

## MARSHALL COMPACTOR, Triple

Code : A017



- Used to determine the resistance of asphalt sample to plastic flow. The test is performed in a fully automatic mechanism and programmed to stop automatically at the required number of blows.
- Supplied with the standard height and rammer weight according to ASTM/EN standards.
- Equipped with a blow counter, wooden pedestal and quick action clamping.

### **STANDARDS**

ASTM D6926 • EN 12697-10 • AASHTO t245

### **TECHNICAL SPECIFICATIONS**

- Consists of 3-mould structure that can compact 3 specimens at the same time
- Blow Rate:  $55 \pm 5$  blows per minute

- Rammer Diameter: 100.33 mm  $\pm$  0.25 mm
- Rammer Weight: 4536 g  $\pm$  9 g
- Drop Height: 457.205 mm  $\pm$  0.635 mm
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

Item	Code
TRIPLE MARSHALL COMPACTOR, AUTOMATIC, STATIONARY BASE - ASTM	A017M04AT
TRIPLE MARSHALL COMPACTOR, AUTOMATIC, STATIONARY BASE - ASTM [60 Hz]	A017M04AK