

ASPHALT > CORING

CORE DRILLING MACHINE (Trailer Mounted)

Code : A065



- Used to take core samples from asphalt.
- Installed on a trailer for easy transportation.
- Equipped with a petrol engine, four stabilizing feet, 100 lt water tank.
- Core Bits for the machine should be ordered separately

STANDARDS

ASTM D 5361 • EN 12697-27

TECHNICAL SPECIFICATIONS

- Trailer Mounted

- Petrol Engine
- Stabilizing Feet (4 ea)
- Water Tank (100 Lt)

ORDERING INFORMATION

| Item | Code |
|---|-----------|
| CORE DRILLING MACHINE - 4 HP (Trailer Mounted) | A065X00XP |
| CORE BIT for ASPHALT PAVEMENT - Ø 50 mm x 400 mm | A065P001H |
| CORE BIT for ASPHALT PAVEMENT - Ø 75 mm x 400 mm | A065P002H |
| CORE BIT for ASPHALT PAVEMENT - Ø 100 mm x 400 mm | A065P003H |
| CORE BIT for ASPHALT PAVEMENT - Ø 150 mm x 400 mm | A065P004H |
| CORE BIT for ASPHALT PAVEMENT - Ø 200 mm x 400 mm | A065P005H |
| FRAME for CORE DRILLING MACHINE (A065) | A065P006H |
| PETROL ENGINE for CORE DRILLING MACHINE | A065P001P |

OTHER PHOTOS

