

# ROLLING STRAIGHTEDGE, TRAVELLING BEAM (High-Low Detector)

Code : A069



- Used for checking planeness/irregularities in concrete and bituminous road surfaces.
- The travelling beam is 3 meters long.
- Deviation of the surface is shown on a scale;
  - in increments of 2 mm from 0 - 10 mm
  - in increments of 5 mm from 10 - 25 mm

## **STANDARDS**

British Department of Transport

## **TECHNICAL SPECIFICATIONS**

- Length : 3 m
- Equipped with scale (up to 25 mm)

## **ORDERING INFORMATION**

Item	Code
ROLLING STRAIGHTEDGE, TRAVELLING BEAM (High-Low Detector) - 3 m	A069X03XH