

CONCRETE > SELF-CONSOLIDATING CONCRETE

V-FUNNEL TEST SET

Code : B025



- The V-Funnel test is a widely used method for assessing the flowability and viscosity of self-compacting concrete (SCC).
- This test helps determine the time it takes for the SCC to flow through a V-shaped funnel, which is indicative of the concrete's ability to pass through tight spaces and around reinforcement without segregation or blockage.

STANDARDS

EN 12350-9

TECHNICAL SPECIFICATIONS

- Made of stainless steel

-
- Stand-mounted

ORDERING INFORMATION

Item	Code
V-FUNNEL TEST SET	B025X00XH