

CONCRETE > WORKABILITY & CONSISTENCY

## K-SLUMP TESTER (Pocket Concrete Flow Meter)

Code : B031



- To assess the level of compaction and the workability of freshly mixed concrete.
- It is utilized for on-site measurements or within testing molds, with test outcomes potentially aligned with slump measurements.

### **STANDARDS**

ASTM C1362

### **ORDERING INFORMATION**

| Item | Code |
|------|------|
|------|------|

K-SLUMP TESTER (Pocket Concrete Flow Meter)

B031X00AH