

CEMENT > WORKABILITY & CONSISTENCY

MINI-SLUMP SET

Code : C020



- The Mini Slump Set is designed to evaluate the workability and consistency of fresh mortar and cement paste, particularly suitable for mixtures with low aggregate content.
- This compact apparatus is ideal for laboratory settings where standard slump tests are impractical due to limited sample volumes or specific testing requirements.
- Constructed from durable materials, the set ensures longevity and consistent performance across repeated testing procedures.
- The mini slump cone facilitates the assessment of flow characteristics, aiding in the optimization of mix designs and quality control processes.
- Its straightforward operation allows for quick and efficient testing, providing immediate insights into the workability of mortar and cementitious mixtures.

TECHNICAL SPECIFICATIONS

- Mini Slump Cone Height: 150 mm
- Top Diameter: 50 mm
- Bottom Diameter: 100 mm
- Tamping Rod: Stainless steel, Ø25 mm, length 190 mm

SET

- Mini slump cone
- Stainless steel tamping rod
- Stainless steel flow table

ORDERING INFORMATION

Item	Code
MINI-SLUMP SET	C020X00XH