

TRICHLOROETHYLENE

Code : G801



- Trichloroethylene (TCE) in engineering testing laboratories is commonly used as a solvent for cleaning equipment, dissolving fats, oils, and greases, and as a degreasing agent in various industrial applications. It can also be used in chemical analysis and testing, although it requires careful handling due to its toxic and carcinogenic properties.

ORDERING INFORMATION

Item	Code
TRICHLOROETHYLENE	G801X00XH
METHYLENE BLUE POWDER - 100 g	G802X00XH
SULPHUR	G803X00XH
GRAPHITE	G804X00XH

SODIUM SULFATE	G805X00XH
----------------	-----------

MAGNESIUM SULFATE	G806X00XH
-------------------	-----------

BARIUM CHLORIDE	G807X00XH
-----------------	-----------

CALCIUM CHLORIDE	G808X00XH
------------------	-----------

SODIUM HYDROXIDE	G809X00XH
------------------	-----------

FORMALDEHYDE (Stock Solution)	G810X00XH
-------------------------------	-----------

SODIUM HEXAMETAPHOSPHATE	G811X00XH
--------------------------	-----------

SILICA GEL	G812X00XH
------------	-----------

PARAFIN WAX	G813X00XH
-------------	-----------

KAOLINITE	G814X00XH
-----------	-----------

QUANTAB CHLORIDE TITRATOR 1175 TYPE	G815X00XH
-------------------------------------	-----------

QUANTAB CHLORIDE TITRATOR 1176 TYPE	G816X00XH
-------------------------------------	-----------