

# CALIBRATION WEIGHT

Code : G950



- Calibration weights are essential tools for the routine verification and adjustment of laboratory balances and scales, ensuring measurement accuracy and compliance with international standards.
- Regular use of these weights facilitates the detection of deviations in balance performance, allowing for timely maintenance and calibration to uphold the integrity of measurement results.
- Sold either separately or as a set of 2110 g, consisting of the following:
  - 1000 g x 1
  - 500 g x 1
  - 200 g x 2
  - 100 g x 1
  - 50 g x 1
  - 20 g x 2
  - 10 g x 1
  - 5 g x 1

- 2 g x 2
- 1 g x 1

## **STANDARDS**

ASTM E617 • ISO 17025

## **ORDERING INFORMATION**

Item	Code
CALIBRATION WEIGHT - 1 g	G950X001H
CALIBRATION WEIGHT - 5 g	G950X005H
CALIBRATION WEIGHT - 10 g	G950X010H
CALIBRATION WEIGHT - 1 kg	G950X01KH
CALIBRATION WEIGHT - 20 g	G950X020H
CALIBRATION WEIGHT - 2 kg	G950X02KH
CALIBRATION WEIGHT - 50 g	G950X050H
CALIBRATION WEIGHT - 5 kg	G950X05KH
CALIBRATION WEIGHT - 100 g	G950X100H
CALIBRATION WEIGHT - 10 kg	G950X10KH
CALIBRATION WEIGHT - 20 kg	G950X20KH
CALIBRATION WEIGHT - 25 kg	G950X25KH
CALIBRATION WEIGHT - 500 g	G950X500H
CALIBRATION WEIGHT - 50 kg	G950X50KH
CALIBRATION WEIGHT SET - 2110 g	G951X001H