


# VIBRATORY COMPACTION HAMMER

Code : T129



- Used to compact concrete, proctor, CBR and bituminious mixtures.
- Considered as an alternative method in determining the relationship between dry density and water content in the compacted soil specimens.

## **STANDARDS**

EN 13286-4 • BS 1377-4 • BS 1924-2

## **SUPPLIED WITH**

- Vibrating Hammer (220 - 240 V / 50 - 60 Hz)
- Supporting frame
- Compaction Foots (Ø102 mm & Ø146 mm)

**ORDERING INFORMATION**

| Item                                                        | Code      |
|-------------------------------------------------------------|-----------|
| VIBRATORY COMPACTION HAMMER (Kango Hammer)                  | T129X00XU |
| TAMPING FOOT - Ø 100 mm for VIBRATORY COMPACTION HAMMER     | T129P001H |
| TAMPING FOOT - Ø 150 mm for VIBRATORY COMPACTION HAMMER     | T129P002H |
| SUPPORTING FRAME for VIBRATORY COMPACTION HAMMER            | T129P003H |
| TAMPING FOOT - 150 x 150 mm for VIBRATORY COMPACTION HAMMER | T129P004H |