

# DESSICATOR

Code : W008



- A desiccator in engineering testing laboratories is used to store samples or materials in a dry environment by absorbing moisture from the air. It helps prevent sample contamination and ensures accurate results in tests that are sensitive to humidity, such as moisture content analysis or when working with hygroscopic materials.

## ORDERING INFORMATION

| Item                       | Code      |
|----------------------------|-----------|
| DESSICATOR - 200 mm        | W008G296H |
| DESSICATOR - 250 mm        | W008G344H |
| DESSICATOR - 300 mm        | W008G420H |
| VACUUM DESSICATOR - 200 mm | W009G296H |

---

VACUUM DESSICATOR - 250 mm

W009G344H

VACUUM DESSICATOR - 300 mm

W009G420H