

HAIR BRUSH

Code : G719



- A hair brush in engineering testing laboratories is used for gently cleaning fine mesh sieves without damaging them. It helps remove fine particles stuck in the mesh, ensuring accurate sieve analysis and maintaining the integrity of the sieve.

ORDERING INFORMATION

Item	Code
WIRE BRUSH	G718X001H
HAIR BRUSH - No: 1	G719X001H
HAIR BRUSH - No: 2	G719X002H
HAIR BRUSH - No: 3	G719X003H
HAIR BRUSH - No: 4	G719X004H

BRUSH with HANDLE

G719X005H