

## SAMPLE TRAY

Code : G768



- A sample tray in engineering testing laboratories is used to hold, transport, and organize material samples during testing or analysis. It helps prevent contamination, ensures proper handling, and allows for easy access and examination of samples.

### ORDERING INFORMATION

| Item                                          | Code      |
|-----------------------------------------------|-----------|
| SAMPLE TRAY - X Small Size, 200 x 300 x 40 mm | G768X001H |
| SAMPLE TRAY - Small Size, 230 x 330 x 40 mm   | G768X002H |
| SAMPLE TRAY - Medium Size, 260 x 380 x 40 mm  | G768X003H |
| SAMPLE TRAY - Large Size, 280 x 430 x 40 mm   | G768X004H |
| SAMPLE TRAY - X Large Size, 300 x 480 x 40 mm | G768X005H |