

# EVAPORATING DISH

Code : W024



- An evaporating dish in engineering testing laboratories is used for evaporating solvents from a solution, leaving behind the solid residue. It is typically made of porcelain or heat-resistant glass, designed to withstand high temperatures during the evaporation process, commonly used in chemical analyses or material testing.

## ORDERING INFORMATION

| Item                        | Code      |
|-----------------------------|-----------|
| EVAPORATING DISH - Ø 30 mm  | W024C015H |
| EVAPORATING DISH - Ø 70 mm  | W024C050H |
| EVAPORATING DISH - Ø 80 mm  | W024C075H |
| EVAPORATING DISH - Ø 100 mm | W025C140H |

EVAPORATING DISH - Ø 125 mm

W025C300H