

ASPHALT > CORING

CORE DRILLING MACHINE (Portable Unit)

Code : A066



- Used to take core samples from asphalt and concrete.
- The machine is designed to be portable and easily carried in a pick-up car.
- Equipped with a petrol engine, four stabilizing screw stands and a vertical screw feed.
- Suitable for Easy Single Person Operation.
- Core Bits for the machine should be ordered separately

STANDARDS

ASTM D 5361 • EN 12697-27

TECHNICAL SPECIFICATIONS

- Portable unit.

- Equipped with:
 - Petrol engine
 - Four stabilizing screws.
 - Vertical screw feed
 - Water inlet
- Engine Configuration: Vertical Shaft
- Engine Fuel: Gasoline
- Starter: Rewind
- Bit Capacity: Up to 6" (150 mm) diameter bits

ORDERING INFORMATION

Item	Code
CORE DRILLING MACHINE (Portable Unit)	A066X00XP
CORE BIT for ASPHALT PAVEMENT - Ø 50 mm x 400 mm	A065P001H
CORE BIT for ASPHALT PAVEMENT - Ø 75 mm x 400 mm	A065P002H
CORE BIT for ASPHALT PAVEMENT - Ø 100 mm x 400 mm	A065P003H
CORE BIT for ASPHALT PAVEMENT - Ø 150 mm x 400 mm	A065P004H
CORE BIT for ASPHALT PAVEMENT - Ø 200 mm x 400 mm	A065P005H
PETROL ENGINE for CORE DRILLING MACHINE	A065P001P

OTHER PHOTOS

