

ASPHALT > BITUMEN PROPERTIES

RING AND BALL APPARATUS

Code : A472



Used to determine the softening point of bitumen immersed in distilled water.

STANDARDS

AASHTO T53 • ASTM D36 • EN 1427

SET

- Hotplate
- Glass Beaker - 600 ml
- Brass frame
- Rings w/ centering guide (2 ea)
- Balls (2 ea)
- Glass Thermometer 200°C

ORDERING INFORMATION

Item	Code
RING and BALL SOFTENING POINT APPARATUS - ASTM	A472H00AU
RING and BALL SOFTENING POINT APPARATUS - EN	A472H00EU
AUTOMATIC RING AND BALL APPARATUS	A472A00XU
STEEL BALL for SOFTENING POINT APPARATUS	A472P001H
BALL CENTERING GUIDE for SOFTENING POINT APPARATUS	A472P002H
BRASS RING for SOFTENING POINT APPARATUS	A472P003H