

## CUBE MOULD

Code : B011



- Cube moulds are essential for preparing concrete specimens for compressive strength testing, ensuring accurate and repeatable results in laboratory and field settings.
- Available in plastic, steel, and cast iron, these moulds cater to various testing needs, offering durability and ease of use depending on the material choice.
- Designed to produce specimens of standard sizes, such as 100 mm, 150 mm, and 200 mm cubes, aligning with international testing standards.
- The robust construction of these moulds ensures they maintain dimensional stability under repeated use, providing consistent specimen shapes for reliable testing outcomes.
- Cast iron mould features a smooth internal surface to facilitate easy demoulding and to produce specimens with accurate dimensions and surface finish.

### **STANDARDS**

ASTM C109 • ASTM C192 • ASTM C617 • ASTM C403 • AASHTO T23 • AASHTO T126 • AASHTO T97 • EN 12390-1 •  
BS 1881:108

## TECHNICAL SPECIFICATIONS

- Available in Plastic, Steel and Cast Iron
- Available Sizes: 100 mm, 150 mm, 200 mm cube dimensions
- Configuration: Single-gang and multi-gang options
- Surface Finish: Smooth interior surfaces for accurate specimen formation

## ORDERING INFORMATION

| Item                                           | Code      |
|------------------------------------------------|-----------|
| CUBE MOULD - Steel (10 cm) : 2-Gang            | B011P001H |
| CUBE MOULD - Cast (15 cm)                      | B011P002H |
| CUBE MOULD - Cast (20 cm)                      | B011P003H |
| CUBE MOULD - Plastic (10 cm) - 2-Gang          | B011P004H |
| CUBE MOULD - Plastic (15 cm)                   | B011P005H |
| COMPACTING BAR (Square Cross-Section) - Steel  | B011P006H |
| COMPACTING BAR (Square Cross-Section) - Wooden | B011P007H |
| PLUG for PLASTIC CUBE MOULD                    | B011P008H |

## OTHER PHOTOS

