

CURING TANK

Code : B140



- Used to cure cylinders, cubes, beams and other samples that requires total immersion in water at a specified constant temperature.
- Equipped with a digital thermostat and indicator.
- Supplied with a base rack and a water circulation pump.

STANDARDS

ASTM C31 • BS 1881-11

TECHNICAL SPECIFICATIONS

- 220 - 240 V / 50 - 60 Hz

SUPPLIED WITH

- Heater Unit with Digital Thermostat and indicator
- Base Rack
- Water Circulation Pump

ORDERING INFORMATION

Item	Code
CURING TANK CIRCULATION PUMP	B140P001U
CURING TANK CONTROLLING UNIT with HEATER	B140P002U
BASE RACK 95 x 195 cm for CURING TANK	B140P003H
CURING TANK - Steel	B140X001U
CURING TANK - Plastic - Large	B140X002U
CURING TANK - Plastic	B140X003U
CURING TANK - Stainless Steel	B140X004U

OTHER PHOTOS

