

SILICA GEL

Code : G812



- Silica gel in engineering testing laboratories is commonly used as a desiccant to absorb moisture and keep samples, equipment, or materials dry. It is also used in chromatography for separating chemical components and as a drying agent in various testing procedures to prevent contamination from moisture.

ORDERING INFORMATION

Item	Code
TRICHLOROETHYLENE	G801X00XH
METHYLENE BLUE POWDER - 100 g	G802X00XH
SULPHUR	G803X00XH
GRAPHITE	G804X00XH
SODIUM SULFATE	G805X00XH

MAGNESIUM SULFATE	G806X00XH
BARIUM CHLORIDE	G807X00XH
CALCIUM CHLORIDE	G808X00XH
SODIUM HYDROXIDE	G809X00XH
FORMALDEHYDE (Stock Solution)	G810X00XH
SODIUM HEXAMETAPHOSPHATE	G811X00XH
SILICA GEL	G812X00XH
PARAFIN WAX	G813X00XH
KAOLINITE	G814X00XH
QUANTAB CHLORIDE TITRATOR 1175 TYPE	G815X00XH
QUANTAB CHLORIDE TITRATOR 1176 TYPE	G816X00XH