

# MOBILE LABORATORY

Code : P001



- We can offer a complete line of Mobile Laboratory with special design to meet all possible environments such as cold or tropical climate.
- The main tree types are:
  - Container
  - Trailer Mounted
  - Van Mounted

## **TECHNICAL SPECIFICATIONS**

- Office / Bedroom section
- Lighting & Electric sockets
- Electric generator
- Air Conditioning / Aspirator Fan

- Water Tank, WC & Sink

## **ORDERING INFORMATION**

| Item                                  | Code      |
|---------------------------------------|-----------|
| MOBILE LABORATORY - 20 DC             | P001X00XX |
| TRAILER for MOBILE LABORATORY - 20 DC | P001P001X |

## **OTHER PHOTOS**

