


MANHOLE COVER TESTER

Code : U007



- The ALFA Manhole Cover Tester is designed for evaluating the load-bearing capacity of manhole covers, gratings, and similar components.
- The system features a high-rigidity steel column frame, ensuring stable and precise testing of samples up to **1200 x 1000 mm** in size. With a maximum capacity of **1000 kN**, the tester is suitable for both routine quality control and advanced research applications.
- The integrated **Touchscreen TFT Control Unit** allows for easy test setup and monitoring, while the included **special software** provides control and data acquisition from a PC.

STANDARDS

AASHTO M306 • EN 124 • DIN 19584 • BSI PAS 26

TECHNICAL SPECIFICATIONS

- Maximum Capacity: 1000 kN
- Maximum Sample Size: 1200 x 1000 mm
- Frame Type: High-rigidity steel column frame
- Hydraulic Power Unit: Suitable for repetitive load cycles
- Control Unit: Touchscreen TFT with PC connectivity
- Data Acquisition: Via PC and built-in touchscreen
- Digital Dial Gauge: 25 mm range, 0.01 mm resolution
- Power Supply: 220 - 240 V / 50 - 60 Hz

EQUIPPED WITH

- Special high-rigidity column frame
- Hydraulic power unit for cyclic loading
- Touchscreen TFT control panel
- Reference measurement bar
- Special control and acquisition software

SUPPLIED WITH

- Reference bar with 25 × 0.01 mm digital dial gauge
- Special control software (PC not included)

ORDERING INFORMATION

Item	Code
------	------

MANHOLE COVER TEST MACHINE

U007A001U

POWER PACK (for Manhole Cover Tester) - 400 Bar

U007P001U