

PORCELAIN MORTAR and RUBBER-HEADED PESTLE

Code : W027



- A porcelain mortar and rubber-headed pestle in engineering testing laboratories are used for grinding and mixing materials without causing contamination or damage. The rubber-headed pestle is gentle on the mortar's surface and is ideal for handling delicate or fragile samples, making it useful for preparing samples in chemical, biological, or material testing.

ORDERING INFORMATION

Item	Code
PORCELAIN MORTAR and PESTLE - Ø 80 mm	W026C080H
PORCELAIN MORTAR and PESTLE - Ø 100 mm	W026C100H
PORCELAIN MORTAR and PESTLE - Ø 130 mm	W026C130H
PORCELAIN MORTAR and PESTLE - Ø 160 mm	W026C160H

PORCELAIN MORTAR and PESTLE - Ø 190 mm	W026C190H
--	-----------

PORCELAIN MORTAR and RUBBER-HEADED PESTLE	W027P001H
---	-----------

GLASS MORTAR and PESTLE - Ø 100 mm	W028G100H
------------------------------------	-----------

GLASS MORTAR and PESTLE - Ø 120 mm	W028G120H
------------------------------------	-----------

GLASS MORTAR and PESTLE - Ø 150 mm	W028G150H
------------------------------------	-----------

GLASS MORTAR and PESTLE - Ø 180 mm	W028G180H
------------------------------------	-----------