


# ACCESSORIES for COMPRESSION TESTING MACHINE

Code : B003P



Some of the accessories and spare parts for compression testers series are listed here. Feel free to contact us for further information and guidance.

## **ORDERING INFORMATION**

Item Name	Item Code
FLEXURAL TESTING MACHINE (Rectangular Frame, 1m x 1m Specimen) - 200 kN	B003A201U
FLEXURAL TESTING MACHINE (Rectangular Frame) - 200 kN	B003A20XU
FLEXURAL TESTING MACHINE (Rectangular Frame) - 300 kN	B003A30XU
FLEXURAL TESTING MACHINE (Rectangular Frame) - 600 kN	B003A60XU
FLEXURE APPARATUS for CONCRETE BEAMS	B003P001H

RECTANGULAR COMPRESSION PLATES - 30x40 cm	B003P003H
SPACER DISK - 40 mm	B003P004H
SPACER DISK - 50 mm	B003P005H
SPACER DISK - 100 mm	B003P006H
SPLITTING TENSILE APPARATUS for Ø 100 mm CYLINDRER SPECIMENS	B003P007H
SPLITTING TENSILE APPARATUS for Ø 150 mm CYLINDRER SPECIMENS	B003P008H
SPLITTING TENSILE APPARATUS for 150 mm CUBES	B003P009H
THERMAL PRINTER for COMPRESSION TESTER (B-001 and B-004)	B003P010H
METAL DOOR for 2000 kN COMPRESSION TESTER	B003P011H
METAL DOOR for 3000 kN COMPRESSION TESTER	B003P012H
SPLITTING TENSILE APPARATUS for 100 mm CUBES	B003P013H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for Ø75x150mm CONCRETE CYLINDERS	B003P020H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for Ø75x150mm CONCRETE CYLINDERS	B003P021H
COMPRESSOMETER/EXTENSOMETER (LPDT) for Ø75x150mm CONCRETE CYLINDERS	B003P022H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for Ø100x200mm CONCRETE CYLINDERS	B003P023H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for Ø100x200mm CONCRETE CYLINDERS	B003P024H
COMPRESSOMETER/EXTENSOMETER (LPDT) for Ø100x200mm CONCRETE CYLINDERS	B003P025H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for Ø150x300mm CONCRETE CYLINDERS	B003P026H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for Ø150x300mm CONCRETE CYLINDERS	B003P027H
COMPRESSOMETER/EXTENSOMETER (LPDT) for Ø150x300mm CONCRETE CYLINDERS	B003P028H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for Ø160x320mm CONCRETE CYLINDERS	B003P029H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for Ø160x320mm CONCRETE CYLINDERS	B003P030H
COMPRESSOMETER/EXTENSOMETER (LPDT) for Ø160x320mm CONCRETE CYLINDERS	B003P031H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for 100x100mm CONCRETE CUBES	B003P032H

COMPRESSOMETER/EXTENSOMETER (DIGITAL) for 100x100mm CONCRETE CUBES	B003P033H
COMPRESSOMETER/EXTENSOMETER (LPDT) for 100x100mm CONCRETE CUBES	B003P034H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for 150x150mm CONCRETE CUBES	B003P035H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for 150x150mm CONCRETE CUBES	B003P036H
COMPRESSOMETER/EXTENSOMETER (LPDT) for 150x150mm CONCRETE CUBES	B003P037H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for 200x200mm CONCRETE CUBES	B003P038H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for 200x200mm CONCRETE CUBES	B003P039H
COMPRESSOMETER/EXTENSOMETER (LPDT) for 200x200mm CONCRETE CUBES	B003P040H
PISTON for COMPRESSION TESTER (2000 kN)	B003P041H
PISTON for COMPRESSION TESTER (3000 kN)	B003P042H
THERMAL PAPER ROLL	B003P043H
STAND / BASE for COMPRESSION TESTER (2000 kN)	B003P044H
BLOCK SLIDING RAIL for 30x40 PLATE for B001	B003P045H
PISTON for COMPRESSION TESTER (1500 kN)	B003P046H
SPACER DISK - 45 mm	B003P047H
ADAPTOR for CURBSTONE TEST	B003P048H
CIRCULAR LOWER COMPRESSION PLATE - Ø 28.5 cm	B003P049H
FLEXURAL RECTANGULAR FRAME - 200 kN	B003P20XH
FLEXURAL RECTANGULAR FRAME - 300 kN	B003P30XH
FLEXURAL RECTANGULAR FRAME - 600 kN	B003P60XH

### OTHER PHOTOS

