

CONCRETE PIPE TESTING MACHINE

Code : B008



- Supplied with a graphical LCD Data Acquisition Control System and an ability to store and report the results.
- The calibration for the equipment can be easily done on 5 pts.
- Pipe length and diameter should be specified at the time of inquiry.

STANDARDS

ASTM C301 • ASTM C497 • EN 1916

TECHNICAL SPECIFICATIONS

Capacity (kN)	400	600	1000
Pipe Diameter (mm)	200 - 1500	200 - 2000	200 - 3700

Pipe Length (mm)	500 - 1500	500 - 2000	1000 - 3000
-------------------------	------------	------------	-------------

ORDERING INFORMATION

Item Name	Item Code
CONCRETE PIPE TESTING MACHINE - 400 kN	B008A04XU
CONCRETE PIPE TESTING MACHINE - 600 kN	B008A06XU
CONCRETE PIPE TESTING MACHINE - 1000 kN	B008A10XU
CONCRETE PIPE TESTING MACHINE - 2000 kN	B008A20XU
V-SHAPED SUPPORT - 1 m	B008P001H
SPARE PISTON (1000 kN) for B008A10XU	B008P002H
LOADING BEAM HINGE for B008A10XU	B008P003H