

## **CAPPING SET**

Code : B015



- Comprises a melting pot and a capping apparatus.
- The melting pot is designed to melt capping compounds and features an adjustable temperature control via a thermoregulator.
- A built-in pilot lamp indicates the heating status.
- The pot is constructed with a double-walled frame for enhanced durability and insulation.

### **STANDARDS**

ASTM C31 • ASTM C192 • EN 12390-3 • AASHTO T23

### **SUPPLIED WITH**

- Capping apparatus (Used for Ø 100 mm and Ø 150 mm cylinders)

**TECHNICAL SPECIFICATIONS**

- Melting pot (4.5 lt):
  - Thermostatically Controlled
  - Insulated double-wall
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

Item Name	Item Code
CAPPING FRAME (for Ø 100 mm and Ø 150 mm)	B015P001H
MELTING POT - 5 lt	B015P002U
MELTING POT - 10 lt	B015P003U
CAPPING FRAME (for Ø 110 mm and Ø 160 mm)	B015P004H
CYLINDER CAPPING SET	B015X00XU

**OTHER PHOTOS**

