

CAPPING PADS & RETAINERS

Code : B016



Unbonded capping offers a safe, precise, and cost-effective alternative to hot sulfur capping compounds for compressive strength testing of concrete specimens.

STANDARDS

ASTM C1231 • AASHTO T22 • AASHTO T851

ORDERING INFORMATION

Item Name	Item Code
CAPPING RETAINER - 3"	B016P001H
CAPPING RETAINER - 4"	B016P002H

CAPPING RETAINER - 6"	B016P003H
CAPPING PAD - Ø 2" / 70 Shore	B016X27XH
CAPPING PAD - Ø 3" / 70 Shore	B016X37XH
CAPPING PAD - Ø 4" / 70 Shore	B016X47XH
CAPPING PAD - Ø 6" / 70 Shore	B016X67XH

OTHER PHOTOS

