

CONCRETE > SELF-CONSOLIDATING CONCRETE

## SLUMP-FLOW TEST SET (SCC)

Code : B021



- The Self-Compacting Concrete (SCC) Slump-Flow Test Set is an essential tool for assessing the flowability and workability of SCC.
- This test determines the ability of the concrete to flow under its own weight without segregation, ensuring that it can fill formwork and encapsulate reinforcement without the need for mechanical vibration.

### STANDARDS

EN 12350-8

### TECHNICAL SPECIFICATIONS

- Set comprises of:
  - G-020/G/C Slump Cone

- Steel Plate (900 x 900 mm) with marked at (Ø 210 mm and Ø 500 mm)

**ORDERING INFORMATION**

Item Name	Item Code
BASE PLATE for SCC SLUMP-FLOW	B021P001H
SCC SLUMP-FLOW TEST SET	B021X00EH