

CONCRETE > SELF-CONSOLIDATING CONCRETE

L-BOX APPARATUS (SCC)

Code : B026



- The L-Box test is a key method for evaluating the passing ability and flow characteristics of self-compacting concrete (SCC).
- This test assesses the concrete's ability to flow through obstacles, such as reinforcement bars, without segregation or blockage, ensuring that it can achieve full compaction under its own weight.

TECHNICAL SPECIFICATIONS

- Made of stainless steel
- Guillotine-form gate

SUPPLIED WITH

- *2 different obstacles:*
 - *2 x Ø12 mm smooth bars having 59 mm gaps*
 - *3 x Ø12 mm smooth bars having 41 mm gaps*

ORDERING INFORMATION

Item Name	Item Code
SLIDING GATE for L-BOX TEST APPARATUS	B026P001H
BAR BARRIER - 3 x Ø 12 mm	B026P002H
BAR BARRIER - 2 x Ø 12 mm	B026P003H
L BOX TEST APPARATUS - Acrylic Glass	B026X001H
L BOX TEST APPARATUS	B026X00XH