

CONCRETE > SELF-CONSOLIDATING CONCRETE

SIEVE SEGREGATION TEST SET (SCC)

Code : B023



- The Sieve Segregation Test Set is utilized to evaluate the segregation resistance of self-compacting concrete (SCC).
- Segregation in SCC can lead to uneven distribution of aggregates, resulting in reduced strength and durability of the final structure. This test helps ensure the homogeneity and quality of the SCC mix.

TECHNICAL SPECIFICATIONS

- Set comprises of:
 - Perforated Sieve - 5 mm Aperture (Ø 300 mm)
 - Pan (Ø 300 mm)
 - Bucket

ORDERING INFORMATION

Item Name	Item Code
SIEVE SEGREGATION TEST SET	B023X00EH