

CONCRETE > SELF-CONSOLIDATING CONCRETE

CONCRETE BLEEDING TEST SET

Code : B028



- The Concrete Bleeding Test Set, compliant with EN 480-4 standards, is essential for assessing the bleeding characteristics of freshly mixed concrete.
- Bleeding refers to the migration of water to the surface of concrete, which can impact its quality, durability, and workability.

STANDARDS

EN 480-4

TECHNICAL SPECIFICATIONS

- Cylindrical container with Cover (Ø 250 x 280 mm)
 - Tamping Rod
 - Pipette (50 ml)
 - Glass Beaker (100 ml)

ORDERING INFORMATION

Item Name	Item Code
CONCRETE BLEEDING TEST SET	B028X00EH