

# VIBRATING TABLE

Code : B125



- Designed from robust steel and has a very steady frame
- Used to compact the concrete specimens in the laboratory
- Equipped with a vibrating motor
- Supplied with a clamp assembly for moulds.

## **STANDARDS**

ASTM C192 • EN 12390-2

## **ORDERING INFORMATION**

Item Name	Item Code
CLAMPING APPARATUS for VIBRATING TABLE (Small)	B125P001H

---

CLAMPING APPARATUS for VIBRATING TABLE (Large)	B125P002H
VIBRATING TABLE (Small)	B125X001U
VIBRATING TABLE (Large)	B125X002U
VIBRATING TABLE (Large) - 380V	B125X002V