

# CURING TANK

Code : B140



- Used to cure cylinders, cubes, beams and other samples that requires total immersion in water at a specified constant temperature.
- Equipped with a digital thermostat and indicator.
- Supplied with a base rack and a water circulation pump.

## **STANDARDS**

ASTM C31 • BS 1881-11

## **TECHNICAL SPECIFICATIONS**

- 220 - 240 V / 50 - 60 Hz

**SUPPLIED WITH**

- Heater Unit with Digital Thermostat and indicator
- Base Rack
- Water Circulation Pump

**ORDERING INFORMATION**

Item Name	Item Code
CURING TANK CIRCULATION PUMP	B140P001U
CURING TANK CONTROLLING UNIT with HEATER	B140P002U
BASE RACK - 90 x 190 cm for CURING TANK	B140P003H
HEATING ELEMENT for CURING TANK (L-Shaped 600 mm x 400 mm) - 1500 W	B140P004H
CURING TANK - Steel (200 x 60 x 63 (h) cm)	B140X001U
CURING TANK - Plastic (200 x 100 x 80 (h) cm)	B140X002U
CURING TANK - Plastic (150 x 100 x 80 (h) cm)	B140X003U
CURING TANK - Stainless Steel (200 x 60 x 63 (h) cm)	B140X004U

**OTHER PHOTOS**

