

CONCRETE > SETTING TIME BY PENETRATION

CONCRETE POCKET PENETROMETER

Code : B304



- Used to determine the initial set of concrete mortar

STANDARDS

ASTM C403 • AASHTO T197 • UNI 7123

TECHNICAL SPECIFICATIONS

- Capacity: 4.5 MPa
- Plunger's surface area: 32.30 mm² (~ 1/20 sq.in)

ORDERING INFORMATION

Item Name	Item Code
CONCRETE POCKET PENETROMETER - Pen-Type	B304X00XH